


Measurements In a Turbulent Channel Flow

Application Note StereoPIV-005


This application describes high spatial resolution 3-component measurements made in a fully developed, turbulent channel flow.

The objective of the current study was to identify, quantify, and understand the dominant structures in the centerplane of a channel flow. A high-resolution StereoPIV system was employed to examine the size, strength and orientation of typical structures in the flow.

The PowerView™ StereoPIV system consisted of a pair of frequency-doubled Nd:YAG lasers, an articulating light delivery arm, light sheet optics, a pair of PIVCAM 10-30 cameras with $1K \times 1K$ pixel resolution, a pair of high-speed interface boards and INSIGHTUltra™ image acquisition and analysis software. The two-camera system using the optimum Scheimpflug configuration was set up in the horizontal plane to make three-component velocity measurements



W-component of velocity



Vorticity distribution

The measurements provided 12, 100 instantaneous three-component velocity vectors per set of frame-straddled, stereoscopic images. Data for 986 statistically independent, three-component velocity fields were captured and computed to compute detailed flow statistics.



TSI Incorporated – 500 Cardigan Road, Shoreview, MN 55126 U.S.A

USA	Tel: +1 800 874 2811	E-mail: fluid@tsi.com	Website: www.tsi.com
UK	Tel: +44 149 4 459200	E-mail: tsiuk@tsi.com	Website: www.tsiinc.co.uk
France	Tel: +33 491 11 87 64	E-mail: tsifrance@tsi.com	Website: www.tsiinc.fr
Germany	Tel: +49 241 523030	E-mail: tsigmbh@tsi.com	Website: www.tsiinc.de
India	Tel: +91 80 41132470	E-mail: tsi-india@tsi.com	
China	Tel: +86 10 8260 1595	E-mail: tsibeijing@tsi.com	
Singapore	Tel: +65 6595 6388	E-mail: tsi-singapore@tsi.com	



TRUST. SCIENCE. INNOVATION.

Contact your local TSI Distributor or visit our website www.tsi.com for more detailed specifications.