

Single-Box Scanning Mobility Particle Sizer™ Spectrometer

Model 3034

Standard replacement parts for Model 3034 Single-Box Scanning Mobility Particle Sizer™ Spectrometer consist of:

<u>Part Number</u>	<u>Part Description</u>
1031580	3034 Maintenance Kit: Includes (2) 1602051; (2) 1602094; (2) 1030013; (2) fan media
1030013	Dacron Screen
1031097	Cyclone (1 L/min; 0.8 Micron)
1036015R	Zero Filter Assembly
1602051	HEPA Capsule Filter (need 2)
1602094	Filter Dispos Inline Ny .6 µm (need 2)
1031486	CPC Fill/Drain Replacement Kit CPC (includes fill bottle, drain bottle, and bottle bracket)

For pricing information, or to place an order you can 1) send the necessary information to answers@tsi.com or 2) contact your local TSI representative (see contact information below) and speak to one of our customer specialists. Please include your contact information (telephone number and mailing address) with all requests.

Scanning Mobility Particle Sizer™ is a trademark of TSI Incorporated.

TSI Incorporated – 500 Cardigan Road, Shoreview, MN 55126 U.S.A

USA Tel: +1 800 874 2811

UK Tel: +44 149 4 459200

France Tel: +33 491 95 21 90

Germany Tel: +49 241 523030

India Tel: +91 80 41132470

China Tel: +86 10 8260 1595

E-mail: particle@tsi.com

E-mail: tsiuk@tsi.com

E-mail: tsifrance@tsi.com

E-mail: tsigmbh@tsi.com

E-mail: tsi-india@tsi.com

E-mail: tsibeijing@tsi.com

Website: www.tsi.com

Website: www.tsiinc.co.uk

Website: www.tsiinc.fr

Website: www.tsiinc.de



TRUST. SCIENCE. INNOVATION.

Contact your local TSI Distributor or visit our website www.tsi.com for more detailed specifications.