

Fluidized Bed Aerosol Generator

Model 3400A

Standard replacement parts for Model 3400A Fluidized Bed Aerosol Generator consist of:

<u>Part Number</u>	<u>Part Description</u>
1031583	3400A Maintenance Kit: Includes (4ft) 1502528; (2) 1502574; (1) Bead Supply Nozzle; (1) Replacement light; (1) Porous platew/sticker; (1) chain pliers; (1) nozzle removal tool; (1) insertion tool
1502528	Bead Chain Gear - bronze
1502563	Replacement Bead Chain, (10-7/16-in long)
1502574	Replacement 100 µm Bronze Beads (comes in 3 lb. jar, 250 ml)
1030658R	Replacement Control Porous Plate w/sticker
1502573	Pliers, Chain

For pricing information, or to place an order you can 1) send the necessary information to answers@tsi.com or 2) contact your local TSI representative (see contact information below) and speak to one of our customer specialists. Please include your contact information (telephone number and mailing address) with all requests.

TSI Incorporated – 500 Cardigan Road, Shoreview, MN 55126 U.S.A

USA Tel: +1 800 874 2811

E-mail: particle@tsi.com

Website: www.tsi.com

UK Tel: +44 149 4 459200

E-mail: tsiuk@tsi.com

Website: www.tsiinc.co.uk

France Tel: +33 491 95 21 90

E-mail: tsifrance@tsi.com

Website: www.tsiinc.fr

Germany Tel: +49 241 523030

E-mail: tsigmbh@tsi.com

Website: www.tsiinc.de

India Tel: +91 80 41132470

E-mail: tsi-india@tsi.com

China Tel: +86 10 8260 1595

E-mail: tsibeijing@tsi.com



TRUST. SCIENCE. INNOVATION.

Contact your local TSI Distributor or visit our website www.tsi.com for more detailed specifications.